


Anti-Mt-Nd4 Rabbit Polyclonal Antibody For Wb, Ihc-P, If/Icc, Elisa - P03905

Item Number: CM0024506



Introduction

Rabbit polyclonal antibody against MT-ND4 (NADH dehydrogenase subunit 4) for WB, IHC-P, IF/ICC, ELISA. Reacts with human, mouse, rat. Suitable for mitochondrial Complex I research.

[Learn More](#)

Parameter	Value
Product Name	MT-ND4 Rabbit pAb
Remarks/Alias	MTND4; ND4; MT-ND4
Species	Human
GeneID	4538
Immunogen	A synthetic peptide corresponding to a sequence within amino acids 401-459 of human MT-ND4 (YP_002124311.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P03905
Protein Weight	52kDa
Shipping	Ice bag