


Anti-Mt-Nd4L Rabbit Polyclonal Antibody For Wb, If/Icc, Elisa - P03901

Item Number: CM0020069



Introduction

Rabbit polyclonal antibody against human MT-ND4L (mitochondrial NADH dehydrogenase subunit 4L). Suitable for WB, IF/ICC, ELISA; cross-reacts with mouse and rat. Useful for investigating Complex I and mitochondrial diseases.

[Learn More](#)

Parameter	Specification
Product Name	MT-ND4L Rabbit pAb
Remarks/Alias	MTND4L; ND4L; MT-ND4L
Species	Human
Gene ID (Human)	4539
Gene ID	4539
Immunogen	Synthetic peptide corresponding to a sequence within amino acids 40 to the C-terminus of human MT-ND4L (AXQ02532.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Mouse, Rat
SWISS	P03901
Protein Weight	11kDa
Shipping	Ice bag