


Anti-Mt-Nd5 Rabbit Polyclonal Antibody For Wb And Elisa - P03915

Item Number: CM0014111



Introduction

Rabbit polyclonal antibody targeting NADH dehydrogenase subunit 5 (MT-ND5) validated for Western blot and ELISA. Cross-reacts with human, mouse, and rat. Immunogen derived from mouse MT-ND5 peptide.

[Learn More](#)

Parameter	Value
Product Name	MT-ND5 Rabbit pAb
Remarks/Alias	ND5; MT-ND5
Species	Mouse
GeneID (Human)	4540
GeneID	17721
Immunogen	A synthetic peptide corresponding to a sequence within amino acids 300-400 of mouse MT-ND5 (NP_904338.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P03915
Protein Weight	67kDa
Shipping	Ice bag