


## Anti-Mt-Nd5 Rabbit Polyclonal Antibody For Wb, Ihc-P, Elisa - P03915

Item Number: CM0021821



### Introduction

Rabbit polyclonal antibody reactive against human MT-ND5, designed for WB, IHC-P, and ELISA applications. Cross-reacts with mouse and rat. Supports research into mitochondrial function, complex I assembly, and mitochondrial diseases.

[Learn More](#)

| Parameter        | Value                                                                                                      |
|------------------|------------------------------------------------------------------------------------------------------------|
| Product Name     | MT-ND5 Rabbit pAb                                                                                          |
| Remarks/Alias    | MTND5; ND5; MT-ND5                                                                                         |
| Species          | Human                                                                                                      |
| GeneID (Human)   | 4540                                                                                                       |
| GeneID           | 4540                                                                                                       |
| Immunogen        | Synthetic peptide corresponding to a sequence within amino acids 500-600 of human MT-ND5 (YP_003024036.1). |
| Source           | Rabbit                                                                                                     |
| Category         | Polyclonal Antibodies                                                                                      |
| Application      | WB, IHC-P, ELISA                                                                                           |
| Cross Reactivity | Human, Mouse, Rat                                                                                          |
| SWISS            | P03915                                                                                                     |
| Protein Weight   | 67kDa                                                                                                      |
| Shipping         | Ice bag                                                                                                    |