


Anti-Mt-Nd5 Rabbit Polyclonal Antibody For Wb, Ihc-P, If-P, Elisa - P03915

Item Number: CM0029054



Introduction

Rabbit polyclonal antibody against MT-ND5 (NADH dehydrogenase subunit 5) for Western blot, immunohistochemistry, immunofluorescence, and ELISA. Cross-reacts with human, mouse, and rat. Ideal for mitochondrial Complex I studies.

[Learn More](#)

Parameter	Value
Product Name	MT-ND5 Rabbit pAb
Remarks/Alias	ND5; MT-ND5
Species	Mouse
Gene ID (Human)	4540
Gene ID (Mouse)	17721
Immunogen	Synthetic peptide A synthetic peptide corresponding to a sequence within amino acids 500-607 of mouse MT-ND5 (NP_904338.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF-P, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P03915
Protein Weight	67kDa
Shipping	Ice bag