


Anti-Mt-Nd6 Rabbit Monoclonal Antibody For Wb, Ihc-P, Elisa - P03923

Item Number: CM0011829



Introduction

High-quality rabbit monoclonal antibody targeting human MT-ND6 (NADH dehydrogenase subunit 6), suitable for WB, IHC-P, and ELISA applications. Detects a 19 kDa protein, essential for mitochondrial Complex I research and disease studies.

[Learn More](#)

Parameter	Value
Product Name	MT-ND6 Rabbit mAb
Remarks/Alias	MTND6; ND6
Species	Human
GeneID (Human)	4541
GeneID	4541
Immunogen	Synthetic peptide corresponding to a sequence within amino acids 134-150 of human MT-ND6 (YP_003024037.1)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human
SWISS	P03923
Protein Weight	19 kDa
Shipping	Ice bag