


## Anti-Mt-Nd6 Rabbit Polyclonal Antibody For Wb, Ihc-P, Elisa - P03923

Item Number: CM0018801



### Introduction

Rabbit polyclonal antibody targeting human MT-ND6, a subunit of mitochondrial complex I. Validated in WB, IHC-P, and ELISA; cross-reacts with mouse and rat. Suitable for mitochondrial disease research.

[Learn More](#)

| Parameter        | Value                                                                                  |
|------------------|----------------------------------------------------------------------------------------|
| Product Name     | MT-ND6 Rabbit pAb                                                                      |
| Remarks/Alias    | MTND6; ND6; MT-ND6                                                                     |
| Species          | Human                                                                                  |
| GeneID (Human)   | 4541                                                                                   |
| GeneID           | 4541                                                                                   |
| Immunogen        | Synthetic peptide corresponding to amino acids 80-174 of human MT-ND6 (YP_003024037.1) |
| Source           | Rabbit                                                                                 |
| Category         | Polyclonal Antibodies                                                                  |
| Application      | WB, IHC-P, ELISA                                                                       |
| Cross Reactivity | Human, Mouse, Rat                                                                      |
| SWISS            | P03923                                                                                 |
| Protein Weight   | 19kDa                                                                                  |
| Shipping         | Ice bag                                                                                |