

Anti-Mt2A Rabbit Polyclonal Antibody For Wb, Ihc-P, Elisa - P02795

Item Number: CM0020420



Introduction

This rabbit polyclonal antibody targets human MT2A (Metallothionein-2A), validated for WB, IHC-P, and ELISA. Cross-reacts with mouse and rat. Ideal for metal homeostasis and detoxification research.

[Learn More](#)

Parameter	Value
Product Name	MT2A Rabbit pAb
Remarks/Alias	MT2, MT-2, MT-II, MT2A
Species	Human
GeneID (Human)	4502
GeneID	4502
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-61 of human MT2A (NP_005944.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P02795
Protein Weight	6kDa
Shipping	Ice bag