


Anti-Mta1 Rabbit Polyclonal Antibody For Wb And Elisa - Q13330

Item Number: CM0012986



Introduction

Rabbit polyclonal antibody against human MTA1, validated for Western blot and ELISA. Detects human and mouse MTA1 (~81 kDa), a NuRD complex component and transcriptional coregulator linked to cancer metastasis. Ideal for cancer biology and gene regulation studies.

[Learn More](#)

Parameter	Value
Product Name	MTA1 Rabbit pAb
Remarks/Alias	MTA1
Species	Human
GeneID	9112
Immunogen	A synthetic peptide corresponding to a sequence within amino acids 600-700 of human MTA1 (NP_004680.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse
SWISS	Q13330
Protein Weight	81kDa
Shipping	Ice bag