

Anti-Mta3 Polyclonal Antibody For Wb And Elisa - Q9Btc8

Item Number: CM0024228

Introduction

Rabbit polyclonal antibody against human MTA3. Validated for WB and ELISA, cross-reacts with human and mouse. Targets a key component of the NuRD complex involved in chromatin remodeling and epithelial architecture.

[Learn More](#)

Parameter	Value
Product Name	MTA3 Rabbit pAb
Remarks/Alias	MTA3
Species	Human
GeneID (Human)	57504
GeneID	57504
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 40-330 of human MTA3 (NP_065795.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse
SWISS	Q9BTC8
Protein Weight	68kDa
Shipping	Ice bag