


Anti-Mtap Polyclonal Antibody For Wb, If/Icc, Elisa - Q13126

Item Number: CM0023523



Introduction

Rabbit polyclonal antibody against MTAP (S-methyl-5'-thioadenosine phosphorylase). Validated for WB, IF/ICC, ELISA. Reacts with human, mouse, rat. Ideal for polyamine metabolism and methionine salvage research.

[Learn More](#)

Parameter	Value
Product Name	MTAP Rabbit pAb
Remarks/Alias	BDMF; MSAP; DMSFH; LGMBF; DMSMFH; c86fus; HEL-249; MTAP
Species	Human
GeneID	4507
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-283 of human MTAP (NP_002442.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q13126
Protein Weight	31kDa
Shipping	Ice bag