


Anti-Mtap Rabbit Polyclonal Antibody For Wb, Ihc-P, If/Icc, Elisa - Q13126

Item Number: CM0017787



Introduction

MTAP rabbit polyclonal antibody validated for WB, IHC-P, IF/ICC, and ELISA. Detects MTAP in human, mouse, and rat. Useful for studying methionine salvage, polyamine metabolism, and cancer.

[Learn More](#)

Parameter	Value
Product Name	MTAP Rabbit pAb
Remarks/Alias	BDMF; MSAP; DMSFH; LGMBF; DMSMFH; c86fus; HEL-249; MTAP
Species	Human
GeneID	4507
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-283 of human MTAP (NP_002442.2).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q13126
Protein Weight	31kDa
Shipping	Ice bag