

Anti-Mtch1 Polyclonal Antibody For Wb, Ihc-P, Elisa - Q9Nzj7

Item Number: CM0028515

Introduction

High-quality rabbit polyclonal antibody against human MTCH1 (mitochondrial carrier homolog 1). Validated for WB, IHC-P, ELISA with cross-reactivity to mouse and rat. Ideal for apoptosis and mitochondrial protein insertion research.

[Learn More](#)

Parameter	Value
Product Name	MTCH1 Rabbit pAb
Remarks/Alias	PSAP; PIG60; CGI-64; SLC25A49; MTCH1
Species	Human
GeneID (Human)	23787
GeneID	23787
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 100-250 of human MTCH1 (NP_001258570.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q9NZJ7
Protein Weight	42kDa
Shipping	Ice bag