

Anti-Mtch2 Polyclonal Antibody For Wb, Ihc-P, If/Icc, Elisa - Q9Y6C9

Item Number: CM0014617

Introduction

Rabbit polyclonal antibody against human MTCH2, a mitochondrial outer membrane insertase and apoptosis receptor. KO validated; suitable for WB, IHC-P, IF/ICC, and ELISA. Cross-reacts with human and mouse samples.

[Learn More](#)

Parameter	Value
Product Name	[KO Validated] MTCH2 Rabbit pAb
Remarks/Alias	MIMP; HSPC032; SLC25A50; H2
Species	Human
GeneID (Human)	23788
GeneID	23788
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 30-175 of human MTCH2 (NP_055157.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, IF/ICC, ELISA
Cross Reactivity	Human, Mouse
SWISS	Q9Y6C9
Protein Weight	33 kDa
Shipping	Ice bag