

Anti-Mtco1 Rabbit Monoclonal Antibody For Wb, Ihc-P, Elisa - P00395

Item Number: CM0007474



Introduction

Rabbit monoclonal antibody targeting human cytochrome c oxidase subunit 1 (MTCO1/COX1). Suitable for Western blot, IHC-paraffin, and ELISA. For mitochondrial respiratory chain research.

[Learn More](#)

Parameter	Value
Product Name	MTCO1 Rabbit mAb
Remarks/Alias	CoxI; COX1; MTCO1
Species	Human
Gene ID (Human)	4512
Gene ID	4512
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 479-513 of human MTCO1 (YP_003024028.1)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human
SWISS	P00395
Protein Weight	57 kDa
Shipping	Ice bag