

Anti-Mtco1 Rabbit Polyclonal Antibody For Wb, Ihc-P, If/Icc, Elisa - P00395

Item Number: CM0019510



Introduction

Rabbit polyclonal antibody against human MTCO1 (Cytochrome c oxidase subunit 1), validated for WB, IHC-P, IF/ICC, ELISA. Suitable for human, mouse, and rat samples. Ideal for mitochondrial disease and bioenergetics research.

[Learn More](#)

Parameter	Value
Product Name	MTCO1 Rabbit pAb
Remarks/Alias	COI; MTCO1; COX1
Species	Human
GeneID	4512
Immunogen	Synthetic peptide corresponding to a sequence within amino acids 400-500 of human MTCO1 (YP_002124304.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P00395
Protein Weight	57kDa
Shipping	Ice bag