

Anti-Mtco2 (Cox2) Rabbit Monoclonal Antibody For Wb - P00403

Item Number: CM0006156



Introduction

Rabbit monoclonal antibody against MTCO2 (COX2), reactive in human, mouse, rat. Validated for WB, IF/ICC, IHC-P, ELISA. Targets mitochondrial cytochrome c oxidase subunit 2, essential for oxidative phosphorylation.

[Learn More](#)

Parameter	Value
Product Name	MTCO2 Rabbit mAb
Remarks/Alias	COX2
Species	Human
GeneID (Human)	17709
GeneID	17709
Immunogen	Synthetic peptide (a synthetic peptide corresponding to a sequence within amino acids 101-200 of mouse MTCO2 (NP_904331.1))
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IF/ICC, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P00403
Protein Weight	26 kDa
Shipping	Ice bag