

## Anti-Mtco2 Polyclonal Antibody For Wb, Ihc-P, If/Icc, Elisa - P00403

Item Number: CM0017767



### Introduction

High-quality rabbit polyclonal antibody against MTCO2 (Cytochrome c oxidase subunit 2) for WB, IHC-P, IF/ICC, ELISA. Validated for mouse and rat, 26 kDa protein. Suitable for mitochondrial function research.

[Learn More](#)

| Parameter        | Value                                                                                                 |
|------------------|-------------------------------------------------------------------------------------------------------|
| Product Name     | MTCO2 Rabbit pAb                                                                                      |
| Remarks/Alias    | COX2; MTCO2                                                                                           |
| Species          | Mouse                                                                                                 |
| GeneID (Human)   | 4513                                                                                                  |
| GeneID           | 17709                                                                                                 |
| Immunogen        | Synthetic peptide corresponding to a sequence within amino acids 101-200 of mouse MTCO2 (NP_904331.1) |
| Source           | Rabbit                                                                                                |
| Category         | Polyclonal Antibodies                                                                                 |
| Application      | WB, IHC-P, IF/ICC, ELISA                                                                              |
| Cross Reactivity | Mouse, Rat                                                                                            |
| SWISS            | P00403                                                                                                |
| Protein Weight   | 26kDa                                                                                                 |
| Shipping         | Ice bag                                                                                               |