

Anti-Mtco2 Rabbit Monoclonal Antibody For Wb, If, Ihc-P, Elisa - P00403

Item Number: CM0011500



Introduction

High-quality rabbit monoclonal antibody targeting human MTCO2 (cytochrome c oxidase subunit 2). Validated for WB, IF/ICC, IHC-P, and ELISA. Suitable for mitochondrial research and oxidative phosphorylation studies.

[Learn More](#)

Parameter	Value
Product Name	MTCO2 Rabbit mAb
Remarks/Alias	COII; MTCO2; COX2
Species	Human
Gene ID	4513
Immunogen	Synthetic Peptide A synthetic peptide corresponding to a sequence within amino acids 128-227 of human MTCO2 (P00403).
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IF/ICC, IF-P, IHC-P, ELISA
Cross Reactivity	Human
SWISS	P00403
Protein Weight	26kDa
Shipping	Ice bag