


Anti-Mtdh Rabbit Polyclonal Antibody For Wb, Ihc-P, If/Icc, Elisa - Q86Ue4

Item Number: CM0028884



Introduction

High-affinity Anti-MTDH (AEG-1/LYRIC) rabbit polyclonal antibody validated for WB, IHC-P, IF/ICC, ELISA; reacts with human, mouse, rat. Essential for cancer metastasis, chemoresistance, and NF-κB signaling research.

[Learn More](#)

Parameter	Value
Product Name	MTDH Rabbit pAb
Remarks/Alias	3D3; AEG1; AEG-1; LYRIC; LYRIC/3D3; MTDH
Species	Human
Gene ID (Human)	92140
Gene ID	92140
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-260 of human MTDH (NP_848927.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q86UE4
Protein Weight	64kDa
Shipping	Ice bag