


Anti-Mterfd1 Rabbit Polyclonal Antibody For Wb, Elisa - Q96E29

Item Number: CM0025086



Introduction

Rabbit polyclonal antibody targeting MTERFD1 (mTERF3), validated for Western blot and ELISA. Cross-reacts with human and mouse. Suitable for mitochondrial transcription, ribosome assembly, and respiratory complex studies.

[Learn More](#)

Parameter	Value
Product Name	MTERFD1 Rabbit pAb
Remarks/Alias	CGI-12; MTERFD1
Species	Human
GeneID (Human)	51001
GeneID	51001
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 69-300 of human MTERFD1 (NP_057026.3)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse
SWISS	Q96E29
Protein Weight	48kDa
Shipping	Ice bag