

Anti-Mterfd3 Rabbit Polyclonal Antibody For Wb, Elisa - Q49Am1

Item Number: CM0027253

Introduction

High-quality rabbit polyclonal antibody against human MTERFD3 (mTERFL), suitable for Western blotting and ELISA, with cross-reactivity to mouse and rat. Ideal for mitochondrial transcription studies.

[Learn More](#)

Parameter	Value
Product Name	MTERFD3 Rabbit pAb
Remarks/Alias	mTERFL; MTERFD3
Species	Human
GeneID (Human)	80298
GeneID	80298
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-250 of human MTERFD3 (NP_079474.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q49AM1
Protein Weight	44kDa
Shipping	Ice bag