

Anti-Mtfp1 Rabbit Polyclonal Antibody For Wb, Elisa - Q9Udx5

Item Number: CM0031084

Introduction

Rabbit polyclonal antibody targeting human MTFP1 (MTP18), a mitochondrial inner membrane protein involved in fission. Validated for WB and ELISA; cross-reacts with human, mouse, and rat. 18kDa.

[Learn More](#)

Parameter	Value
Product Name	MTFP1 Rabbit pAb
Remarks/Alias	MTP18; HSPC242; MTFP1
Species	Human
GeneID (human)	51537
GeneID	51537
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-166 of human MTFP1 (NP_057582.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q9UDX5
Protein Weight	18kDa
Shipping	Ice bag