


Anti-Mth1 Polyclonal Antibody For Wb, Elisa - P36639

Item Number: CM0030692



Introduction

Rabbit polyclonal antibody against human MTH1 (NUDT1), validated for WB and ELISA. Suitable for exploring DNA damage repair, oxidative stress, and nucleotide pool sanitization pathways. SWISS: P36639.

[Learn More](#)

Parameter	Value
Product Name	MTH1 Rabbit pAb
Remarks/Alias	MTH1
Species	Human
GeneID (Human)	4521
GeneID	4521
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-179 of human MTH1 (NP_945192.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human
SWISS	P36639
Protein Weight	18kDa
Shipping	Ice bag