

Anti-Mthfd1L Polyclonal Antibody For Wb, Ihc-P, Elisa - Q6Ub35

Item Number: CM0030115

Introduction

Monofunctional C1-tetrahydrofolate synthase (MTHFD1L) rabbit polyclonal antibody validated for WB, IHC-P, ELISA in human, mouse, rat. Targets mitochondrial one-carbon metabolism enzyme; useful for cancer and metabolic research.

[Learn More](#)

Parameter	Value
Product Name	MTHFD1L Rabbit pAb
Remarks/Alias	FTHFSDC1; MTC1THFS; dj292B18.2; MTHFD1L
Species	Human
GeneID (human)	25902
GeneID	25902
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 759-978 of human MTHFD1L (NP_056255.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q6UB35
Protein Weight	106kDa
Shipping	Ice bag