

Anti-Mthfd2 Rabbit Monoclonal Antibody For Wb, Elisa - P13995

Item Number: CM0007035



Introduction

This high-quality MTHFD2 rabbit monoclonal antibody validated for WB and ELISA, reacting with human and mouse MTHFD2. A key reagent for mitochondrial one-carbon metabolism studies. (P13995)

[Learn More](#)

Parameter	Details
Product Name	MTHFD2 Rabbit mAb
Remarks/Alias	NMDMC; MTHFD2
Species	Human
GeneID (Human)	10797
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 40-250 of human MTHFD2 (NP_006627.2)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse
SWISS	P13995
Protein Weight	38kDa
Shipping	Ice bag