

Anti-Mthfd2 Rabbit Polyclonal Antibody For Wb, Ihc-P, If/Icc, Elisa - P13995

Item Number: CM0017762



Introduction

Rabbit polyclonal antibody targeting human MTHFD2 (bifunctional methylenetetrahydrofolate dehydrogenase/cyclohydrolase) for WB, IHC-P, IF/ICC, ELISA. Reacts with human, mouse, rat. Ideal for studying mitochondrial one-carbon metabolism.

[Learn More](#)

Parameter	Value
Product Name	MTHFD2 Rabbit pAb
Remarks/Alias	NMDMC; MTHFD2
Species	Human
GeneID (Human)	10797
GeneID	10797
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 40-250 of human MTHFD2 (NP_006627.2).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P13995
Protein Weight	38kDa
Shipping	Ice bag