

Anti-Mthfr Monoclonal Antibody For Wb, If/Icc, Elisa - P42898

Item Number: CM0006677



Introduction

High-quality rabbit monoclonal antibody targeting human MTHFR (methylenetetrahydrofolate reductase). Validated for WB, IF/ICC, ELISA with cross-reactivity to mouse and rat. Suitable for folate pathway research and diagnostic development.

[Learn More](#)

Parameter	Value
Product Name	MTHFR Rabbit mAb
Remarks/Alias	MTHFR
Species	Human
GeneID (Human)	4524
GeneID	4524
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 363-656 of human MTHFR (NP_005948.3)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P42898
Protein Weight	75 kDa
Shipping	Ice bag