


## Anti-Mtif3 Polyclonal Antibody For Wb, Elisa - Q9H2K0

Item Number: CM0029735



### Introduction

MTIF3 rabbit polyclonal antibody validated for WB and ELISA. Specific to human and rat MTIF3, a mitochondrial translation initiation factor that promotes 28S ribosomal subunit availability for protein synthesis.

[Learn More](#)

| Parameter        | Value                                                                                                             |
|------------------|-------------------------------------------------------------------------------------------------------------------|
| Product Name     | MTIF3 Rabbit pAb                                                                                                  |
| Remarks/Alias    | IF3mt; MTIF3                                                                                                      |
| Species          | Human                                                                                                             |
| GeneID (Human)   | 219402                                                                                                            |
| GeneID           | 219402                                                                                                            |
| Immunogen        | Recombinant fusion protein containing a sequence corresponding to amino acids 1-278 of human MTIF3 (NP_690876.3). |
| Source           | Rabbit                                                                                                            |
| Category         | Polyclonal Antibodies                                                                                             |
| Application      | WB, ELISA                                                                                                         |
| Cross Reactivity | Human, Rat                                                                                                        |
| SWISS            | Q9H2K0                                                                                                            |
| Protein Weight   | 32kDa                                                                                                             |
| Shipping         | Ice bag                                                                                                           |