


Anti-Mtif3 Polyclonal Antibody For Wb, Elisa - Q9H2K0

Item Number: CM0016065



Introduction

Rabbit polyclonal antibody against human MTIF3 (IF3mt), validated for WB and ELISA applications. Detects MTIF3 in human, mouse, and rat samples. Ideal for studying mitochondrial protein synthesis.

[Learn More](#)

Parameter	Value
Product Name	MTIF3 Rabbit pAb
Remarks/Alias	IF3mt; MTIF3
Species	Human
Gene ID (Human)	219402
Gene ID	219402
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-278 of human MTIF3 (NP_690876.3).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q9H2K0
Protein Weight	32kDa
Shipping	Ice bag