

Anti-Mtmr1 Rabbit Polyclonal Antibody For Wb And Elisa - Q13613

Item Number: CM0018817



Introduction

Rabbit polyclonal antibody against MTMR1, a phosphatidylinositol-3-phosphate phosphatase, for WB and ELISA in human and mouse, supporting lipid metabolism research.

[Learn More](#)

Parameter	Specification
Product Name	MTMR1 Rabbit pAb
Remarks/Alias	MTMR1
Species	Human
Gene ID (Human)	8776
Gene ID	8776
Immunogen	Synthetic peptide corresponding to amino acids 300-400 of human MTMR1 (NP_003819.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse
SWISS	Q13613
Protein Weight	75kDa
Shipping	Ice bag