

Anti-Mtmr14 Polyclonal Antibody For Wb, If/Icc, Elisa - Q8Nce2

Item Number: CM0029610

Introduction

Rabbit polyclonal antibody against MTMR14, a phosphatidylinositol-3,5-bisphosphate 3-phosphatase that dephosphorylates PI3P and PI(3,5)P2. Validated for WB, IF/ICC, and ELISA; reacts with human, mouse, and rat samples. Suitable for lipid signaling and autophagy research.

[Learn More](#)

Parameter	Value
Product Name	MTMR14 Rabbit pAb
Remarks/Alias	C3orf29; MTMR14
Species	Human
GeneID (human)	64419
GeneID	64419
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 335-500 of human MTMR14 (NP_001070993.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q8NCE2
Protein Weight	72kDa
Shipping	Ice bag