

Anti-Mtmr14 Rabbit Polyclonal Antibody For Wb, Elisa - Q8Nce2

Item Number: CM0021520

Introduction

MTMR14 Rabbit Polyclonal Antibody validated for WB and ELISA, with reactivity against human and mouse MTMR14. This antibody targets the lipid phosphatase MTMR14, involved in phosphatidylinositol metabolism, making it suitable for research in cell signaling and disorders linked to phosphoinositide regulation.

[Learn More](#)

Parameter	Value
Product Name	MTMR14 Rabbit pAb
Remarks/Alias	C3orf29; MTMR14
Species	Human
Gene ID (Human)	64419
Gene ID	64419
Immunogen	Recombinant protein Recombinant fusion protein containing a sequence corresponding to amino acids 1-280 of human MTMR14 (NP_001070993.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse
SWISS	Q8NCE2
Protein Weight	72kDa
Shipping	Ice bag