


Anti-Mtmr2 Rabbit Polyclonal Antibody For Wb, Elisa - Q13614

Item Number: CM0025431



Introduction

Rabbit polyclonal antibody against human MTMR2, validated for WB and ELISA. Cross-reacts with mouse and rat. Immunogen corresponds to amino acids 404-643. Ideal for autophagy and myotubularin-related protein research.

[Learn More](#)

Parameter	Value
Product Name	MTMR2 Rabbit pAb
Remarks/Alias	CMT4B; CMT4B1; MTMR2
Species	Human
GeneID	8898
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 404-643 of human MTMR2 (NP_057240.3)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q13614
Protein Weight	73kDa
Shipping	Ice bag