

Anti-Mtmr4 Rabbit Polyclonal Antibody For Wb, If/Icc, Elisa - Q9Nya4

Item Number: CM0016918

Introduction

High-quality rabbit polyclonal antibody against human MTMR4, a phosphatidylinositol-3,5-bisphosphate 3-phosphatase involved in endosomal trafficking and signaling. Validated for WB, IF/ICC, ELISA with cross-reactivity to mouse and rat.

[Learn More](#)

Parameter	Value
Product Name	MTMR4 Rabbit pAb
Remarks/Alias	ZFYVE11; FYVE-DSP2; MTMR4
Species	Human
GeneID (Human)	9110
GeneID	9110
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 920-1195 of human MTMR4 (NP_004678.3)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q9NYA4
Protein Weight	133kDa
Shipping	Ice bag