


Anti-Mtmr6 Rabbit Polyclonal Antibody For Wb, Elisa - Q9Y217

Item Number: CM0015981



Introduction

Rabbit polyclonal antibody targeting MTMR6, a lipid phosphatase involved in T-cell signaling and endocytosis. Validated for WB and ELISA. Reacts with human, mouse, and rat. Suitable for immunodetection and research.

[Learn More](#)

Parameter	Value
Product Name	MTMR6 Rabbit pAb
Remarks/Alias	MTMR6
Species	Human
Gene ID (Human)	9107
Gene ID	9107
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 382-621 of human MTMR6 (NP_004676.3)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q9Y217
Protein Weight	72kDa
Shipping	Ice bag