

Anti-Mtmr9 Rabbit Polyclonal Antibody For Wb, Ihc-P, If/Icc, Elisa - Q96Qg7

Item Number: CM0016675

Introduction

Rabbit polyclonal antibody against MTMR9, validated for WB, IHC-P, IF/ICC, ELISA in human, mouse, rat. Ideal for studying myotubularin-related phosphatase regulation, macroautophagy, and apoptosis. - Q96QG7

[Learn More](#)

Parameter	Value
Product Name	MTMR9 Rabbit pAb
Remarks/Alias	MTMR8; C8orf9; LIP-STYX; MTMR9
Species	Human
Gene ID	66036
Immunogen	Recombinant protein: Recombinant fusion protein containing a sequence corresponding to amino acids 1-250 of human MTMR9 (NP_056273.2).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q96QG7
Protein Weight	63kDa
Shipping	Ice bag