


Anti-Mtnr1A Rabbit Polyclonal Antibody For Wb And Elisa - P48039

Item Number: CM0027676



Introduction

MTNR1A Rabbit Polyclonal Antibody (P48039) validated for western blot (WB) and ELISA. Reacts with human and rat melatonin receptor type 1A. Ideal for circadian rhythm, sleep, and melatonin receptor studies.

[Learn More](#)

Parameter	Value
Product Name	MTNR1A Rabbit pAb
Remarks/Alias	MT1; MEL-1A-R
Species	Human
GeneID (Human)	4543
GeneID	4543
Immunogen	Synthetic peptide corresponding to a sequence within amino acids 261-350 of human MTNR1A (NP_005949.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Rat
SWISS	P48039
Protein Weight	39kDa
Shipping	Ice bag