


## Anti-Mtnr1A Rabbit Polyclonal Antibody (Wb, Ihc-P, Elisa) - P48039

Item Number: CM0013748



### Introduction

Rabbit polyclonal antibody against human MTNR1A (melatonin receptor type 1A), validated for WB, IHC-P, ELISA. Cross-reacts with mouse and rat. Predicted MW 39kDa. Ideal for circadian rhythm and sleep research.

[Learn More](#)

| Parameter        | Value                                                                                                               |
|------------------|---------------------------------------------------------------------------------------------------------------------|
| Product Name     | MTNR1A Rabbit pAb                                                                                                   |
| Remarks/Alias    | MT1; MEL-1A-R; MTNR1A                                                                                               |
| Species          | Human                                                                                                               |
| Gene ID          | 4543 (Human)                                                                                                        |
| Immunogen        | Recombinant fusion protein containing a sequence corresponding to amino acids 281-350 of human MTNR1A (NP_005949.1) |
| Source           | Rabbit                                                                                                              |
| Category         | Polyclonal Antibodies                                                                                               |
| Application      | WB, IHC-P, ELISA                                                                                                    |
| Cross Reactivity | Human, Mouse, Rat                                                                                                   |
| SWISS            | P48039                                                                                                              |
| Protein Weight   | 39kDa                                                                                                               |
| Shipping         | Ice bag                                                                                                             |