

Anti-Mtor Rabbit Polyclonal Antibody For Wb, If/Icc, Ip, Elisa - P42345

Item Number: CM0014884



Introduction

Rabbit polyclonal antibody targeting human and mouse mTOR (P42345). Offered for Western blot, immunofluorescence/immunocytochemistry, immunoprecipitation, and ELISA. Ideal for research on mTOR signaling in cell growth, metabolism, and autophagy.

[Learn More](#)

Parameter	Value
Product Name	mTOR Rabbit pAb
Remarks/Alias	SKS; FRAP; FRAP1; FRAP2; RAFT1; RAPT1; mTOR
Species	Human
Gene ID	2475
Immunogen	Synthetic peptide corresponding to a sequence within amino acids 2400-2500 of human mTOR (NP_004949.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, IP, ELISA
Cross Reactivity	Human, Mouse
SWISS	P42345
Protein Weight	289kDa
Shipping	Ice bag