

Anti-Mtpap Rabbit Polyclonal Antibody For Wb, Elisa - Q9Nvv4

Item Number: CM0019844



Introduction

Rabbit polyclonal antibody against human MTPAP (PAPD1/SPAX4/TENT6), validated for WB and ELISA. Reacts with human and mouse. Suitable for mitochondrial RNA processing and histone mRNA degradation studies.

[Learn More](#)

Parameter	Value
Product Name	MTPAP Rabbit pAb
Remarks/Alias	PAPD1; SPAX4; TENT6; MTPAP
Species	Human
Gene ID (Human)	55149
Gene ID	55149
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 383-582 of human MTPAP (NP_060579.3).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse
SWISS	Q9NVV4
Protein Weight	66kDa
Shipping	Ice bag