


Anti-Mtrr Rabbit Polyclonal Antibody For Wb, Elisa - Q9UbK8

Item Number: CM0013798



Introduction

Rabbit polyclonal antibody against human MTRR (methionine synthase reductase), validated for Western blot (WB) and ELISA. Suitable for studying methionine and folate homeostasis. Reacts with human samples; predicted molecular weight 78kDa.

[Learn More](#)

Parameter	Specification
Product Name	MTRR Rabbit pAb
Remarks/Alias	MSR; cbIE; MTRR
Species	Human
GeneID (Human)	4552
GeneID	4552
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 50-230 of human MTRR (NP_002445.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human
SWISS	Q9UBK8
Protein Weight	78kDa
Shipping	Ice bag