


Anti-Mtss1 Rabbit Polyclonal Antibody For Wb And Elisa - O43312

Item Number: CM0017350



Introduction

High-quality rabbit polyclonal antibody against human MTSS1 (MIM, metastasis suppressor 1) validated for Western blot and ELISA applications. Reacts with human, mouse, and rat MTSS1. Ideal for cancer metastasis and cytoskeleton research.

[Learn More](#)

Parameter	Value
Product Name	MTSS1 Rabbit pAb
Remarks/Alias	MIM; MIMA; MIMB; MTSS1
Species	Human
GeneID	9788
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 630-730 of human MTSS1 (NP_055566.3)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	O43312
Protein Weight	82kDa
Shipping	Ice bag