

Anti-Mttp Polyclonal Antibody For Wb, Elisa - P55157

Item Number: CM0017676



Introduction

Anti-MTTP polyclonal antibody (pAb) validated for WB and ELISA in human, mouse, rat. Targets microsomal triglyceride transfer protein large subunit, key for lipoprotein assembly. Ideal for metabolism studies.

[Learn More](#)

Parameter	Value
Product Name	MTTP Rabbit pAb
Remarks/Alias	ABL; MTP; MTTP
Species	Human
Gene ID	4547
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 19-160 of human MTTP (NP_000244.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P55157
Protein Weight	99kDa
Shipping	Ice bag