

Anti-Mttp Rabbit Polyclonal Antibody For Wb, If/Icc, Elisa - P55157

Item Number: CM0024520



Introduction

Rabbit polyclonal antibody targeting human MTP (microsomal triglyceride transfer protein large subunit). Validated for WB, IF/ICC, ELISA; cross-reacts with mouse, rat. Ideal for research on lipid transport and lipoprotein assembly.

[Learn More](#)

Parameter	Value
Product Name	MTP Rabbit pAb
Remarks/Alias	ABL; MTP; MTP
Species	Human
GeneID (Human)	4547
GeneID	4547
Immunogen	Synthetic peptide corresponding to a sequence within amino acids 550-650 of human MTP (NP_000244.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Mouse, Rat
SWISS	P55157
Protein Weight	99kDa
Shipping	Ice bag