

## Anti-Mtus1 Polyclonal Antibody For Wb, Elisa - Q9Uld2

Item Number: CM0015590

### Introduction

Rabbit polyclonal antibody targeting human MTUS1, validated for Western blot and ELISA. MTUS1 is a microtubule-associated tumor suppressor involved in cell proliferation and mitosis; suitable for oncology research.

[Learn More](#)

| Parameter        | Value                                                                                                              |
|------------------|--------------------------------------------------------------------------------------------------------------------|
| Product Name     | MTUS1 Rabbit pAb                                                                                                   |
| Remarks/Alias    | ATBP, ATIP, ICIS, MP44, ATIP3, MTSG1, MTUS1                                                                        |
| Species          | Human                                                                                                              |
| Gene ID (Human)  | 57509                                                                                                              |
| Gene ID          | 57509                                                                                                              |
| Immunogen        | Recombinant fusion protein containing a sequence corresponding to amino acids 207-436 of human MTUS1 (NP_065800.1) |
| Source           | Rabbit                                                                                                             |
| Category         | Polyclonal Antibodies                                                                                              |
| Application      | WB, ELISA                                                                                                          |
| Cross Reactivity | Human                                                                                                              |
| SWISS            | Q9ULD2                                                                                                             |
| Protein Weight   | 141kDa                                                                                                             |
| Shipping         | Ice bag                                                                                                            |