

Anti-Mtus2 Polyclonal Antibody (Tip150) For Wb & Elisa - Q5Jr59

Item Number: CM0020223

Introduction

Rabbit polyclonal antibody targeting human MTUS2 (microtubule-associated tumor suppressor candidate 2, TIP150). Validated for Western blot and ELISA; cross-reacts with mouse and rat. Ideal for microtubule dynamics studies.

[Learn More](#)

Parameter	Value
Product Name	MTUS2 Rabbit pAb
Remarks/Alias	ICIS; CAZIP; TIP150; KIAA0774; MTUS2
Species	Human
GeneID (human)	23281
GeneID	23281
Immunogen	Recombinant Protein Recombinant fusion protein containing a sequence corresponding to amino acids 49-348 of human MTUS2 (NP_056048.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Mouse, Rat
SWISS	Q5JR59
Protein Weight	150kDa
Shipping	Ice bag