

Anti-Mtx2 Polyclonal Antibody For Wb, If/Icc, Elisa - 075431

Item Number: CM0028973

Introduction

CamelBio supplies anti-MTX2 rabbit polyclonal antibody validated for WB, IF/ICC, and ELISA. Detects human, mouse, and rat metaxin-2, a mitochondrial outer membrane protein involved in protein import. Ideal for IVD research and mitochondrial studies.

[Learn More](#)

Parameter	Value
Product Name	MTX2 Rabbit pAb
Remarks/Alias	MDPS; metaxin-2; MTX2
Species	Human
Gene ID (Human)	10651
Gene ID	10651
Immunogen	Recombinant protein (Recombinant fusion protein containing a sequence corresponding to amino acids 1-263 of human MTX2 (NP_006545.1))
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	O75431
Protein Weight	30kDa
Shipping	Ice bag