


Anti-Mu Opioid Receptor (Mor) Rabbit Polyclonal Antibody For Wb, If-P, Ihc-P, Elisa - P35372

Item Number: CM0022588



Introduction

Rabbit polyclonal antibody targeting human Mu Opioid Receptor (MOR). Suitable for WB, IF-P, IHC-P, and ELISA. Reacts with human, mouse, and rat MOR. Key tool for opioid signaling and addiction research.

[Learn More](#)

Parameter	Value
Product Name	Mu Opioid Receptor(MOR) Rabbit pAb
Remarks/Alias	MOP; MOR; LMOR; MOR1; OPRM; M-OR-1; Mu Opioid Receptor(MOR)
Species	Human
Gene ID (Human)	4988
Gene ID	4988
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 354-388 of human Mu Opioid Receptor(MOR) (NP_000905.3)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF-P, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P35372
Protein Weight	45kDa
Shipping	Ice bag