

Anti-Muc1-C Monoclonal Antibody For Wb, Ihc-P, Elisa - P15941

Item Number: CM0003841



Introduction

High-specificity rabbit monoclonal antibody targeting the C-terminal domain of human MUC1 (Mucin-1). Validated for WB, IHC-P, ELISA; cross-reacts with human, mouse, rat. Ideal for cancer biomarker research, cell signaling, and epithelial cell studies.

[Learn More](#)

Parameter	Value
Product Name	MUC1-C Rabbit mAb
Remarks/Alias	EMA; MCD; PEM; PUM; KL-6; MAM6; MCKD; PEMT; CD227; H23AG; MCKD1; MUC-1; ADMCKD; ADTKD2; Ca15-3; ADMCKD1; CA 15-3; MUC-1/X; MUC1/ZD; MUC-1/SEC; MUC1
Species	Human
GenelD (Human)	4582
GenelD	4582
Immunogen	A synthetic peptide corresponding to a sequence within amino acids 1156-1255 of human MUC1 (P15941)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P15941
Protein Weight	122kDa
Shipping	Ice bag